


# Work Package 3: Proton and ion capture

Poram Ruksasakchai  
(On behalf of WP3)

# Gabor lens





Using the electron plasma to produce a lens. It reduces the focal length to 1% of that produced by a magnetic field alone.





Using the Penning-Malmberg trap, a linear array of electrically biased cylinders arranged along the axis of a uniform magnetic field.

*Aymar, G., Becker, T., Boogert, S., Borghesi, M., Bingham, R., Brenner, C., ... & Xiao, R. (2020). LhARA: the laser-hybrid accelerator for radiobiological applications. Frontiers in Physics, 8, 567738.*



# Current Apparatus (preliminary activity)





# Preliminary results (October 2023)





Load





Trap




Hold



Eject & Detect




# Vary the load time




# Temperature of the trapped $e^-$

Measuring the temperature by detecting the number of electrons escaping from the trap when reducing the trap potential.  
Load time: 100ms





Temperature  $\sim 1\text{eV}$  at equilibrium

# Vary hold time



# Rotating Electric fields / Rotating wall



A six-segment rotating wall electrode is used to control plasma radii. The relative phase of the signal applied to each sector of the electrode is labelled.


This gives a rotating electric field perpendicular to the axis of symmetry of the plasma.

↪ The electric field induces an electric dipole moment in the plasma, leading to plasma compression.

<https://alpha.web.cern.ch/science/rotating-wall>




# Rotating wall




Load time: 200ms

RW frequency: 1MHz, amplitude: 0.33V



Hold+RW time: 250ms-10s



# RW + Cooling of electrons with CF<sub>4</sub>




RW frequency: 4MHz, RW amplitude: 0.33V, Vary Hold+RW time: 17us-30.38s




# Next step ...

- Replacing MCP with Phosphor screen.



- Replacing E-gun linear filament with a spiral filament to obtain more symmetric plasma.



- Optimizing experimental parameters to obtain higher number of the electron with longer lifetime.

# WP3 Personnel



W. Bertsche

M. Charlton

S. Eriksson

D. P. van der Werf

C. A. Isaac

C. Baker

P. Ruksasakchai

T. Dascalu

J. Fajans

J. Wurtele

B. Bingham

R. Hugtenburg

J. Purden

A. Knoll

E. Bennet

PDRA, etc.

