

Proton and ion capture

Gabor lens

- The focal length (f) of the Gabor lens:

$$\frac{1}{f} = \frac{e^2 n_e l}{4\epsilon_0 U}$$

where e is electric charge of the electron



n_e is the plasma density

l is the length of the plasma

ϵ_0 is the permittivity of free space


U is the kinetic energy of the positively charge particle.

- Penning-Malmberg trap







Aymar, G., Becker, T., Boogert, S., Borghesi, M., Bingham, R., Brenner, C., ... & Xiao, R. (2020). LhARA: the laser-hybrid accelerator for radiobiological applications. *Frontiers in Physics*, 8, 567738.

Swansea Penning-Malmberg trap





Electron source




Trapping e^- plasmas

- Rotating wall + Cooling gas (CO_2)





September 2024:
 Trap potential: -400V
 Peak density: $1.1 \times 10^{15} \text{ m}^{-3}$
 Radius: $\sim 5 \text{ mm}$










Ions in the trap


Imaging ions in the trap






Measure current of the trapped ions



Eliminate the ions






Peak density





Peak density = $1 \times 10^{15} m^{-3}$

Lifetime of the plasma

(in the raised trap)



Reproducibility and stability



Conclusion and outlook



Density

- The peak density of the plasma is $1 \times 10^{15} \text{ m}^{-3}$.

Radius



- $\sim 5 \text{ mm}$.

Lifetime


- As long as the RW is applied, the plasma remains stable.

Next steps ...

- Investigate and improve the plasma profile.



WP3 Personnel



W. Bertsche
M. Charlton
S. Eriksson
D. P. van der Werf
C. A. Isaac
C. Baker
P. Ruksasakchai

T. Dascalu
J. Fajans
J. Wurtele
B. Bingham
R. Hugtenburg
J. Purden
A. Knoll
E. Bennet

