

Proton and ion capture

Gabor lens

- The focal length (f) of the Gabor lens:

$$\frac{1}{f} = \frac{e^2 n_e l}{4\epsilon_0 U}$$

where e is electric charge of the electron


n_e is the plasma density

l is the length of the plasma


ϵ_0 is the permittivity of free space

U is the kinetic energy of the positively charge particle.



- Penning-Malmberg trap




Aymar, G., Becker, T., Boogert, S., Borghesi, M., Bingham, R., Brenner, C., ... & Xiao, R. (2020). LhARA: the laser-hybrid accelerator for radiobiological applications. *Frontiers in Physics*, 8, 567738.



Experimental setup







Previous results (-140V trap)



To generate a stable and high-density plasma, two key components are required:

- 1. Ionization:** Oscillating electric field (Osc) + Background gas
- 2. Dimple that traps ions**



- Cause dissociation (breaking apart the molecule), generating secondary species like CO, O₂, or carbon-based ions that may affect plasma behaviour and diagnostics.

Documentation of plasma results with/without CO₂ gas



Seed Plasma (without CO₂)

1. N(t): Number of particles vs. load time
2. Lifetime: Number of particles vs. hold time
3. Temperature
4. Mode

Lens Plasma

Cooling gas: CO₂ (fixed pressure)

Oscillating field: applied to E2


1. N(t): Number of particles vs. load time
2. Stability: Number of particles vs. hold time (oscillating field on)
3. T(t): Temperature vs. time during oscillating field application
4. Lifetime: Number of particles vs. hold time (oscillating field off)
5. Mode

Documentation of plasma results with/without CO₂ gas

Seed Plasma (without CO₂)


1. N(t): Number of particles vs. load time
2. Lifetime: Number of particles vs. hold time
3. Temperature
4. Mode

Swansea positron beamline





E-gun: Investigation and improvements

Before:




After:




Findings:

1. Loose Wires
2. Unstable Mounting



Improvements Implemented:

1. Wire Management
2. Enhanced Mechanical Stability







Plasma results without CO₂ gas

Seed Plasma (without CO₂)

1. $N(t)$: Number of particles vs. load time
2. Lifetime: Number of particles vs. hold time
3. Temperature
4. Mode

Experimental sequence (Trap depth: -140V)

- 1) Load: 3 μ s-10s
- 2) Hold: 20 μ s
- 3) Dump + Trigger digitizer/CAM: 5ms







Plasma results without CO₂ gas

Seed Plasma (without CO₂)

1. $N(t)$: Number of particles vs. load time
2. Lifetime: Number of particles vs. hold time
3. Temperature
4. Mode

Experimental sequence (Trap depth: -140V)

- 1) Load: 1s
- 2) Hold: 3 μ s-5s
- 3) Dump + Trigger digitizer/CAM: 5ms







Plasma results without CO₂ gas

Seed Plasma (without CO₂)

1. N(t): Number of particles vs. load time
2. Lifetime: Number of particles vs. hold time
3. Temperature
4. Mode

Experimental sequence (Trap depth: -140V)

- 1) Load: 1s
- 2) Hold: 20 μ s
- 3) Dump + Trigger digitizer/CAM: 5ms



Plasma results with CO₂ gas



Lens Plasma





Cooling gas: CO₂ (fixed pressure)

Oscillating field: applied to E2

1. **N(t): Number of particles vs. load time**
2. Stability: Number of particles vs. hold time (oscillating field on)
3. T(t): Temperature vs. time during oscillating field application
4. Lifetime: Number of particles vs. hold time (oscillating field off)
5. Mode

Experimental sequence (Trap depth: -140V)

- 1) **Load: 3us-10s**
- 2) Hold: 20us
- 3) Dump + Trigger digitizer/CAM: 5ms



Plasma results with CO₂ gas






Lens Plasma

Cooling gas: CO₂ (fixed pressure)

Oscillating field: applied to E2

1. N(t): Number of particles vs. load time
2. Stability: Number of particles vs. hold time (oscillating field on)
3. T(t): Temperature vs. time during oscillating field application
4. Lifetime: Number of particles vs. hold time (oscillating field off)
5. Mode

Plasma results with CO₂ gas



Plasma results with CO₂ gas






Lens Plasma

Cooling gas: CO₂ (fixed pressure)

Oscillating field: applied to E2

1. N(t): Number of particles vs. load time
2. Stability: Number of particles vs. hold time (oscillating field on)
3. T(t): Temperature vs. time during oscillating field application
4. Lifetime: Number of particles vs. hold time (oscillating field off)
5. Mode

Plasma results with CO₂ gas



Conclusion and outlook






To generate a stable and high-density plasma, two key components are required:

1. **Ionization:** Oscillating electric field (Osc) + Background gas
2. **Dimple that traps ions**

Notably, this plasma can be initiated using a very small number of seed electrons.

Next steps ...

- Investigate plasma modes.
- Resume experiments using the high-voltage (HV) amplifiers.
- Conduct the experiments using Ar instead of CO₂.
- Use numerical simulation to verify our diagnostics.



WP3 Personnel



W. Bertsche

M. Charlton

S. Eriksson

D. P. van der Werf

C. A. Isaac

C. Baker

P. Ruksasakchai

T. Dascalu

J. Fajans

J. Wurtele

B. Bingham

R. Hugtenburg

J. Purden

A. Knoll

E. Bennet



SUPERCOMPUTING WALES
UWCHGYFRIFIADURA CYMRU

